

Variety of test Cartridges

	iBG10-OXI*	iBG10*	iBG8-OXI	iBG8	iBC4	iBG3
pH	—	—	—	—		—
pCO ₂	—	—	—	—		—
pO ₂	—	—	—	—		—
Na ⁺	—	—	—	—	—	
K ⁺	—	—	—	—	—	
Cl ⁻	—	—	—	—	—	
Ca ⁺⁺	—	—	—	—	—	
Lac*	—	—				
Glu*	—	—				
tHb	—	—	—	—		—
sO ₂	—		—			
FO ₂ Hb	—		—			
FHHb	—		—			
FCOHb	—		—			
FMetHb	—		—			
Calculated parameters: cH ⁺ , pH(std), HCO ₃ ⁻ (act), HCO ₃ ⁻ (std), tCO ₂ , BE(ecf), BE(act), BE(B), BB(B), ctO ₂ , P50, sO ₂ (est), Ca ⁺⁺ (7.4), AnGap, Hct(est), pO ₂ (A-a), pO ₂ (a/A), RI, pO ₂ /FIO ₂ , cH+(T), pH(T), pCO ₂ (T), pO ₂ (T), pO ₂ (A-a)(T), pO ₂ (a/A)(T), RI(T), pO ₂ (T)/FIO ₂ , mOsm						

*coming soon

Specification

Test Principle	Optical Fluorescence
Throughput	Results in 60s after sample aspiration
Sample volume	65uL / 95uL
Battery	6800mAh
Display	10.1 inch color touch screen, 1280*800
Dimensions(L*W*H)	245*271.4*292.18 mm
Weight	< 5 Kg (battery included)
Operation Environment	10°C-30°C; %RH: 25%-80%; 57 -106.6 KPa

About Edan

Edan is a healthcare company dedicated to improving the human condition around the world by delivering value-driven, innovative and high-quality medical products and services. For over 20 years, Edan has been pioneering a comprehensive line of medical solutions that address a broad range of healthcare practices including:

- Diagnostic ECG
 - Patient Monitoring
 - OB/GYN
- Ultrasound Imaging
 - Point-of-Care Testing
- *In-Vitro* Diagnostics
 - Veterinary

Healthcare professionals around the world depend on Edan's breakthrough medical technologies and outstanding customer support.



ENG-POCT-i20
V1.0-20220711

Global Headquarters:
Edan Instruments, Inc. | 15 Jinhui Road, Pingshan District, Shenzhen 518122 P.R. China | +86.755.26898326 | www.edan.com | info@edan.com

U.S. and Canada inquiries:
EDAN Diagnostics, Inc. | 9918 Via Pasar, San Diego, CA 92126 | +1.858.750.3066 | www.edandiagnostics.com | edan-info@edandiagnostics.com

© Edan Instruments, Inc. All rights reserved. Features and specifications are subject to change without prior notice. No reproduction, copy or transmission may be made without written permission. Not all products or features are available in all countries, contact Edan for local availability.




A Brighter Solution of
Critical Care
i20 Blood Gas and
Chemistry Analysis System



i20 Blood Gas and Chemistry Analysis System

In critical care, blood gases test is often required to be done rapidly with small sample volume to ensure that results are accurate. With years' professional experience in this area, Edan's next-generation - i20, meets all the demands of healthcare providers for blood gases analysis, makes it an ideal blood gases analysis product in optimal patient care and satisfaction!



10.1 inches High-resolution Color touch screen with user-friendly software minimize training requirements

Auto-sampling reduces pre-analytic sample errors

Integrated barcode scanner/thermal printer speeds data management

Build-in battery allows 24-hours standby time

Simple



- Auto-sampling minimizes hands on time and training requirement
- Zero maintenance, portable, ready for use at any time
- Single-use cartridge avoids contamination
- 10.1 inches' color touch screen with user-friendly software

Smart

- Calibrator with build-in NFC chip allows on-board directly
- Short TAT achieves earlier diagnosis and interventions
- Integrated data transmits report to LIS/HIS/DMS synchronically
- Up to 24 hours-standby time supports better flexibility

Superior

- Best in Class accuracy and precision
- Comprehensive test menu, includes BGEM and immunoassay*
- Smaller sample volume, minimum as 65ul
- Larger capacity, up to 60,000 patient data storage



Next-Generation Test Cartridge

- Microfluidics combines with Optical Fluorescence achieve accurate and reliable test result
- More than 40 parameters in one cartridge provides effective critical care to the patients
- Innovative optical fluorescence microchip allows long storage life at room temperature



i20 Quality Control System

Triple Guarantee the Result

